

<u>Vendor Daemon Compatibility When Installing 64-bit</u> <u>ADMET Predictor</u>

Versions Affected

If you previously installed any 32-bit application from Simulations Plus Inc. (GastroPlus, DDDPlus, MembranePlus, MedChemStudio, or PKPlus), and are now installing the 64-bit version of ADMET Predictor, you will need to ensure that the vendor daemon, which handles the licensing communication, is upgraded to the 64-bit version.

You may disregard this notice if you have only installed the 64-bit versions of ADMET Predictor or MedChem Studio, and have never installed 32-bit versions of Simulations Plus software.

Ensuring A Compatible Vendor Daemon

The Simulations Plus vendor daemon is an executable file (simplus.exe) that resides within the Flexera License Server installation directory.

This directory will either be on your local PC system (for Standalone-type installations) or a directory on the file system of your remote license server (for Network-type installations).

The default installation location for the Flexera License Server is here:

C:\Program Files (x86)\FlexNet Publisher License Server Manager\

You may verify the version of the daemon by right clicking on simplus.exe, and choosing "Properties":

		2			
w folder				Open	
Name	Date modified	Туре	0	Run as administrator	
📔 conf	4/8/2016 12:14 PM	File fol	Ŭ	Troubleshoot compatibility	
\mu demo	4/8/2016 12:14 PM	File fol	2	Edit with Notepad++	
퉬 eventlog	4/8/2016 12:14 PM	File fol	۲	Scan with ESET Endpoint Antivirus	
퉬 examples	4/8/2016 12:14 PM	File fol		Advanced options	
🐌 licenses	1/17/2017 10:55 AM	File fol		1	
퉬 uninstall	4/8/2016 12:26 PM	File fol	C	IDrive	•
🐌 web	4/8/2016 12:14 PM	File fol		Add to archive	
퉬 wsdl	4/8/2016 12:14 PM	File fol	۱	Add to "simplus.rar"	
hasp_rt.exe	8/8/2014 3:15 AM	Applic	۱	Compress and email	
🚳 haspsrm_win32.dll	8/8/2014 3:15 AM	Applic		Compress to "simplus.rar" and email	
💷 Imadmin.exe	9/1/2014 4:02 PM	Applic		Pin to Taskbar	
README-swidtag-Imadmin.txt	8/8/2014 3:16 AM	Text D		Pin to Start Menu	
regid.2009-06.com.flexerasoftware_6d08d	8/8/2014 3:16 AM	SWIDT		Restore previous versions	
💷 simplus.exe	6/25/2013 5:12 PM	Applic		Condition 1	
				Send to	1
				PGP Desktop	•
				Cut	
				Сору	
۲. III.				Create shortcut	
/2013 5:12 PM Date created: 8/17/2016 2:22 PM				Delete	
MB				Delete	
			9	Kename	
			>	Properties	

Then, click on the "Details" tab of the properties box:

💷 simplus.exe Proj	simplus.exe Properties			
General Compatib	General Compatibility Security Details Previous Versions			
Property	Value			
Description	Description			
File description	Licensing Daemon x64			
Type File version Product name	Application 1.0.0.0			
Product version Copyright Size	1,0,0,0 Copyright Simulations Plus Inc. All Rights			
Date modified	6/25/2013 5:12 PM			
Remove Propertie	Remove Properties and Personal Information			
	OK Cancel Apply			

If the file description states "**Licensing Daemon x64**," then you already have the 64-bit vendor daemon installed.

If the file description *does not* state "**Licensing Daemon x64**," then you will need to replace your current daemon with the x64 (64-bit) version by following these steps:

- 1. Acquire the x64 licensing daemon.
 - a. You may obtain the x64 licensing daemon from the ADMET Predictor x64 installation package, within the "Vendor Daemon" folder.
 - b. Alternatively, the x64 daemon may be downloaded from the following location:

https://simplus-downloads.com/Utility/daemon_64/simplus.exe

- 2. Copy the 64-bit simplus.exe into the Flexera License Server installation directory, overwriting the current simplus.exe that resides there.
 - a. Your IT department may need to perform this task if you have a Networktype installation with a central license server that only IT personnel have access to.
- 3. Restart the licensing service to initialize the newly-replaced daemon.
 - a. Click on the Start button, type "services.msc" and press Enter.
 - b. The list of installed services should appear. Locate the Imadmin service.
 - c. Right-click on Imadmin and choose Restart.

If you need any further assistance, please feel free to call us at **(888) 266-9294** or email us at <u>support@simulations-plus.com</u>