




SCIENCE + SOFTWARE = SUCCESS

Installation Instructions for ADMET Predictor 13.0

Before you begin

In order to install ADMET Predictor, you must have administrative privileges.


Installation




The first step in the process is to extract the files from the AP13.0_Install.zip file.

Right-click on the AP13.0_Install.zip and select “Extract All...”.

Accept the default destination folder and click on “Extract”.




In the AP13.0_Install directory, right-click on “setup.exe” and select “Run as administrator”.




ADMET Predictor requires several items which may or may not already be installed on your system. These items will be installed as needed.

Click on “Install” and the required items will be installed.




The ADMET Predictor installation wizard will begin.

Click on “Next” to begin.



The license agreement will be displayed.


Select “Yes” if you agree to the licensing terms and continue the installation.



Select the installation directory for ADMET Predictor.


Unless you have a specific need to install in a different directory, the default location is recommended.

Click "Next".




The installation has the information it needs to continue.

Click "Install" to install ADMET Predictor on your computer.




A progress bar will be displayed while ADMET Predictor installs the required files and components.



When the installation is complete, the Finish dialog will be displayed.

Click on “Finish” to complete the installation.

ADMET Predictor is now installed on your computer.



The first time ADMET Predictor starts it will need to know where to find a license file. Here are two options:

1. Provide a path to a Reprise license file, folder or server, or
 2. Provide an activation key.
- Refer to your email from our Licensing Team.

Then click “OK” and ADMET Predictor will run.